

兔抗 GTF2F1 多克隆抗体

中文名称: 兔抗 GTF2F1 多克隆抗体

英文名称: Anti-GTF2F1 rabbit polyclonal antibody

别名: general transcription factor IIF subunit 1; BTF4; RAP74; TF2F1; TFIIF

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: GTF2F1

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	TFIIF is a general transcription initiation factor that binds to RNA polymerase II and helps to recruit it to the initiation complex in collaboration with TFIIB. It promotes transcription elongation.
Applications:	ELISA, IHC
Name of antibody:	GTF2F1
Immunogen:	Fusion protein of human GTF2F1
Full name:	general transcription factor IIF subunit 1
Synonyms:	BTF4; RAP74; TF2F1; TFIIF
SwissProt:	P35269
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer

IHC Recommend dilution:

30-150

