

兔抗 HOMER1 多克隆抗体

- 中文名称：兔抗 HOMER1 多克隆抗体
- 英文名称：Anti-HOMER1 rabbit polyclonal antibody
- 别名：HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 抗原：HOMER1
- 宿主：Rabbit
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the homer family of dendritic proteins. Members of this family regulate group 1 metabotropic glutamate receptor function.
Applications:	ELISA, WB, IHC
Name of antibody:	HOMER1
Immunogen:	Fusion protein of human HOMER1
Full name:	homer homolog 1
Synonyms :	HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C
SwissProt:	Q86YM7

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human lung cancer and human thyroid cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	40 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	500-2000



