

### 兔抗 GTF3C3 多克隆抗体

中文名称：兔抗 GTF3C3 多克隆抗体

英文名称：Anti-GTF3C3 rabbit polyclonal antibody

别名：general transcription factor IIIC subunit 3; TFIIIC102; TFIIICgamma; TFiiC2-102

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GTF3C3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is part of the TFIIIC2 complex, which binds to the promoters of small nuclear and cytoplasmic RNA genes in order to recruit RNA polymerase III. The TFIIIC2 complex is composed of six subunits. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GTF3C3
<b>Immunogen:</b>	Synthetic peptide of human GTF3C3
<b>Full name:</b>	general transcription factor IIIC subunit 3
<b>Synonyms:</b>	TFIIIC102; TFIIICgamma; TFiiC2-102
<b>SwissProt:</b>	Q9Y5Q9

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	40-200



