

兔抗 HOXA1B1D1 多克隆抗体

中文名称：兔抗 HOXA1/B1/D1 多克隆抗体

英文名称：Anti-HOXA1/B1/D1 rabbit polyclonal antibody

别名：BSAS; HOX1; HOX1F/HOX2; HCFP3; HOX2I; Hox-2.9/HOX4; HOX4G; Hox-4.7

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：HOXA1/B1/D1

宿主：Rabbit

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Acts on the anterior body structures.
Applications:	IHC
Name of antibody:	HOXA1/B1/D1
Immunogen:	Synthesized peptide derived from internal of human HOXA1/B1/D1.
Full name:	homeobox A1/B1/D1
Synonyms :	BSAS; HOX1; HOX1F/HOX2; HCFP3; HOX2I; Hox-2.9/HOX4; HOX4G; Hox-4.7

SwissProt:	P14653/P49639/Q9GZZ0
IHC positive control:	Human testis tissue
IHC Recommend dilution:	50-100

