

兔抗 MPV17 多克隆抗体

中文名称：兔抗 MPV17 多克隆抗体

英文名称： Anti-MPV17 rabbit polyclonal antibody

别名： MPV17, mitochondrial inner membrane protein; SYM1; MTDPS6

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MPV17

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a mitochondrial inner membrane protein that is implicated in the metabolism of reactive oxygen species. Mutations in this gene have been associated with the hepatocerebral form of mitochondrial DNA depletion syndrome (MDDS).
Applications:	ELISA, IHC
Name of antibody:	MPV17
Immunogen:	Fusion protein of human MPV17
Full name:	MPV17, mitochondrial inner membrane protein
Synonyms:	SYM1; MTDPS6
SwissProt:	P39210

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	40-200

