

兔抗 GUCA1A 多克隆抗体

中文名称：兔抗 GUCA1A 多克隆抗体

英文名称： Anti-GUCA1A rabbit polyclonal antibody

别名： guanylate cyclase activator 1A; COD3; GCAP; GUCA; GCAP1; GUCA1; CORD14; C6orf131

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GUCA1A

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an enzyme that plays a role in the recovery of retinal photoreceptors from photobleaching. This enzyme promotes the activity of retinal guanylyl cyclase-1 (GC1) at low calcium concentrations and inhibits GC1 at high calcium concentrations. Mutations in this gene can cause cone dystrophy 3 and rod dystrophy 14. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	GUCA1A

Immunogen:	Fusion protein of human GUCA1A
Full name:	guanylate cyclase activator 1A
Synonyms:	COD3; GCAP; GUCA; GCAP1; GUCA1; CORD14; C6orf131
SwissProt:	P43080
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human tonsil
IHC Recommend dilution:	100-300



