

## 兔抗 HOXA3 多克隆抗体

中文名称：兔抗 HOXA3 多克隆抗体

英文名称：Anti-HOXA3 rabbit polyclonal antibody

别名：HOX1; HOX1E

相关类别：一抗

抗原：HOXA3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. Three transcript variants encoding two different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	HOXA3

<b>Immunogen:</b>	Synthetic peptide of human HOXA3
<b>Full name:</b>	homeobox A3
<b>Synonyms :</b>	HOX1; HOX1E
<b>SwissProt:</b>	O43365
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	25-100

