

兔抗 GUCA1B 多克隆抗体

- 中文名称：兔抗 GUCA1B 多克隆抗体
- 英文名称：Anti-GUCA1B rabbit polyclonal antibody
- 别名：guanylate cyclase activator 1B; RP48; GCAP2; GUCA2
- 抗原：GUCA1B
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a calcium-binding protein that activates photoreceptor guanylate cyclases. This gene may have arisen due to a gene duplication event since there is a highly similar gene clustered with it on chromosome 6. Mutations in this gene can cause a form of retinitis pigmentosa.
Applications:	ELISA, WB, IHC
Name of antibody:	GUCA1B
Immunogen:	Synthetic peptide of human GUCA1B
Full name:	guanylate cyclase activator 1B
Synonyms:	RP48; GCAP2; GUCA2

SwissProt:	Q9UMX6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	23 kDa
WB Positive control:	Mouse eye tissue lysate
WB Recommended dilution:	500-2000



