

## 兔抗 HOXA4 多克隆抗体

中文名称：兔抗 HOXA4 多克隆抗体

英文名称：Anti-HOXA4 rabbit polyclonal antibody

别名：HOX1, HOX1D

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HOXA4

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	HOXA4
<b>Immunogen:</b>	Synthetic peptide of human HOXA4
<b>Full name:</b>	homeobox A4

<b>Synonyms :</b>	HOX1, HOX1D
<b>SwissProt:</b>	Q00056
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human liver cancer and Human colon cancer
<b>IHC Recommend dilution :</b>	25-100

