

## 兔抗 GAP43 多克隆抗体

中文名称：兔抗 GAP43 多克隆抗体

英文名称： Anti-GAP43 rabbit polyclonal antibody

别名： B-50; PP46

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GAP43

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene has been termed a 'growth' or 'plasticity' protein because it is expressed at high levels in neuronal growth cones during development and axonal regeneration. This protein is considered a crucial component of an effective regenerative response in the nervous system. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GAP43
<b>Immunogen:</b>	Synthetic peptide of human GAP43

<b>Full name:</b>	growth associated protein 43
<b>Synonyms:</b>	B-50; PP46
<b>SwissProt:</b>	P17677
<b>ELISA Recommended dilution:</b>	500-1000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	25-100

