

兔抗 GUCY1B2 多克隆抗体

中文名称：兔抗 GUCY1B2 多克隆抗体

英文名称： Anti-GUCY1B2 rabbit polyclonal antibody

别名： guanylate cyclase 1 soluble subunit beta 2 (pseudogene); GC-SB2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GUCY1B2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Guanylate Cyclase is an heterodimer of an alpha and a beta chain. Its enzymatic activity transforms GTP to GMP releasing a diphosphate. It binds 1 or 2 hemes per heterodimer. It is activated by nitric oxide in the presence of magnesium or manganese ions. There are two types of guanylate cyclase: soluble form and membrane-associated receptor form.
Applications:	ELISA, IHC
Name of antibody:	GUCY1B2
Immunogen:	Synthetic peptide of human GUCY1B2
Full name:	guanylate cyclase 1 soluble subunit beta 2 (pseudo

	gene)
Synonyms:	GC-SB2
SwissProt:	O75343
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human tonsil
IHC Recommend dilution:	20-100



