

兔抗 HOXA11 多克隆抗体

中文名称：兔抗 HOXA11 多克隆抗体

英文名称：Anti-HOXA11 rabbit polyclonal antibody

别名：HOX1, HOX1I

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HOXA11

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. This gene is involved in the regulation of uterine development and is required for female fertility. Mutations in this gene can cause radio-ulnar synostosis with amegakaryocytic thrombocytopenia.
Applications:	ELISA, WB

Name of antibody:	HOXA11
Immunogen:	Synthetic peptide of human HOXA11
Full name:	homeobox A11
Synonyms :	HOX1, HOX1I
SwissProt:	P31270
ELISA Recommended dilution:	1000-5000
WB Predicted band size:	34 kDa
WB Positive control:	Mouse muscle tissue
WB Recommended dilution:	500-2000

