

HSF1 抗原（重组蛋白）

中文名称：HSF1 抗原（重组蛋白）

英文名称：HSF1 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 371-529 amino acids of human HSF1

技术规格

Full name:	heat shock transcription factor 1
Synonyms:	HSTF1
Swissprot:	Q00613
Gene Accession:	BC014638
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The product of this gene is a heat-shock transcription factor. Transcription of heat-shock genes is rapidly induced after temperature stress. Hsp90, by itself and/or associated with multichaperone complexes, is a major repressor of this gene.