

## 兔抗 HOXB2 多克隆抗体

中文名称：兔抗 HOXB2 多克隆抗体

英文名称：Anti-HOXB2 rabbit polyclonal antibody

别名：K8; HOX2; HOX2H; Hox-2.8

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HOXB2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the Antp homeobox family and encodes a nuclear protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox B genes located on chromosome 17. The encoded protein functions as a sequence-specific transcription factor that is involved in development. Increased expression of this gene is associated with pancreatic cancer.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	HOXB2
<b>Immunogen:</b>	Synthetic peptide of human HOXB2
<b>Full name:</b>	homeobox B2
<b>Synonyms :</b>	K8; HOX2; HOX2H; Hox-2.8

<b>SwissProt:</b>	P14652
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	38 kDa
<b>WB Positive control:</b>	Human bladder cancer tissue
<b>WB Recommended dilution:</b>	200-1000

