

兔抗 GARNL3 多克隆抗体

- 中文名称: 兔抗 GARNL3 多克隆抗体
- 英文名称: Anti-GARNL3 rabbit polyclonal antibody
- 别名: GTPase activating Rap/RanGAP domain like 3; bA356B19.1
- 相关类别: 一抗
- 抗原: GARNL3
- 储存: 冷冻 (-20℃)
- 宿主: Rabbit
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

| | |
|------------------------------------|--|
| Background: | The function of this protein remains unknown. |
| Applications: | ELISA, IHC |
| Name of antibody: | GARNL3 |
| Immunogen: | Synthetic peptide of human GARNL3 |
| Full name: | GTPase activating Rap/RanGAP domain like 3 |
| Synonyms: | bA356B19.1 |
| SwissProt: | Q5VWW2 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human colorectal cancer and Human thyroid cancer |
| IHC Recommend dilution: | 40-200 |



