

IFNG 抗原（重组蛋白）

中文名称：IFNG 抗原（重组蛋白）

英文名称：IFNG Antigen (Recombinant Protein)

别名：IFG; IFI

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 24-161 amino acids of human IFNG

技术规格

Full name:	interferon, gamma
Synonyms:	IFG; IFI
Swissprot:	P01579
Gene Accession:	BC070256
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the type II interferon family. The protein encoded is a soluble cytokine with antiviral, immunoregulatory and anti-tumor properties and is a potent activator of macrophages. Mutations in this gene are associated with aplastic anemia.