

## 兔抗 GAS2L1 多克隆抗体

中文名称：兔抗 GAS2L1 多克隆抗体

英文名称： Anti-GAS2L1 rabbit polyclonal antibody

别名： GAR22

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GAS2L1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the growth arrest-specific 2 protein family. This protein binds components of the cytoskeleton and may be involved in mediating interactions between microtubules and microfilaments. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 9.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GAS2L1
<b>Immunogen:</b>	Synthetic peptide of human GAS2L1
<b>Full name:</b>	growth arrest specific 2 like 1
<b>Synonyms:</b>	GAR22

<b>SwissProt:</b>	Q99501
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	10-50

