

兔抗 GZMB 多克隆抗体

中文名称：兔抗 GZMB 多克隆抗体

英文名称： Anti-GZMB rabbit polyclonal antibody

别名： HLP; CCPI; CGL1; CSPB; SECT; CGL-1; CSP-B; CTLA1; CTSL1

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： GZMB

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Cytolytic T lymphocytes (CTL) and natural killer (NK) cells share the remarkable ability to recognize, bind, and lyse specific target cells. They are thought to protect their host by lysing cells bearing on their surface 'nonself' antigens, usually peptides or proteins resulting from infection by intracellular pathogens. The protein encoded by this gene is crucial for the rapid induction of target cell apoptosis by CTL in cell-mediated immune response.

Applications:	ELISA, IHC
Name of antibody:	GZMB
Immunogen:	Fusion protein of human GZMB
Full name:	granzyme B (granzyme 2, cytotoxic T-lymphocyte-associated serine esterase 1)
Synonyms :	HLP; CCPI; CGL1; CSPB; SECT; CGL-1; CSP-B; CTLA 1; CTSGL1
SwissProt:	P10144
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

