

兔抗 GAS2 多克隆抗体

中文名称：兔抗 GAS2 多克隆抗体

英文名称： Anti-GAS2 rabbit polyclonal antibody

储 存： 冷冻（-20℃）

抗 原： GAS2

宿 主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a caspase-3 substrate that plays a role in regulating microfilament and cell shape changes during apoptosis. It can also modulate cell susceptibility to p53-dependent apoptosis by inhibiting calpain activity. Multiple alternatively spliced variants, encoding the same protein, have been identified.
Applications:	ELISA, WB, IHC
Name of antibody:	GAS2
Immunogen:	Full length fusion protein
Full name:	growth arrest-specific 2
SwissProt:	O43903
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and human gastric cancer

IHC Recommend dilution:	25-100
WB Predicted band size:	35 kDa
WB Positive control:	Human fetal liver tissue
WB Recommended dilution:	200-1000



