

## 兔抗 GAS6 多克隆抗体

中文名称：兔抗 GAS6 多克隆抗体

英文名称：Anti-GAS6 rabbit polyclonal antibody

别名：AXSF, AXLLG

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GAS6

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

This gene product is a gamma-carboxyglutamic acid (Gla)-containing protein thought to be involved in the stimulation of cell proliferation, and may play a role in thrombosis. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. Ligand for tyrosine-protein kinase receptors AXL, TYRO3 and MER whose signaling is implicated in cell growth and survival, cell adhesion and cell migration. GAS6/AXL signaling plays a role in various processes such as endothelial cell survival during acidification by preventing apoptosis, optimal cytokine signaling during human natural killer cell development, hep

	atic regeneration, gonadotropin-releasing hormone neuron survival and migration, platelet activation, or regulation of thrombotic responses.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GAS6
<b>Immunogen:</b>	Synthetic peptide of human GAS6
<b>Full name:</b>	growth arrest-specific 6
<b>Synonyms :</b>	AXSF, AXLLG
<b>SwissProt:</b>	Q14393
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human ovarian cancer and Human cervical cancer
<b>IHC Recommend dilution:</b>	10-50

