

兔抗 GASK1A 多克隆抗体

中文名称：兔抗 GASK1A 多克隆抗体

英文名称：Anti-GASK1A rabbit polyclonal antibody

别名：golgi associated kinase 1A; C3orf41; FAM198A

抗原：GASK1A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	GASK1A
Immunogen:	Synthetic peptide of human GASK1A
Full name:	golgi associated kinase 1A
Synonyms:	C3orf41; FAM198A
SwissProt:	Q9UFP1
IHC positive control:	Human brain
IHC Recommend dilution:	20-100

