

## 兔抗 H1-3 多克隆抗体

- 中文名称： 兔抗 H1-3 多克隆抗体
- 英文名称： Anti-H1-3 rabbit polyclonal antibody
- 别名： H1.3 linker histone, cluster member; H1D; H1.3; H1
- 相关类别： 一抗
- 抗原： H1-3
- 储存： 冷冻（-20℃）
- 宿主： Rabbit
- 反应种属： Human, Mouse
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

### 技术规格

**Background:**

Histones are basic nuclear proteins responsible for nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone

	e H1 family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6. [provided by RefSeq, Aug 2015]
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	H1-3
<b>Immunogen:</b>	Synthetic peptide of human H1-3
<b>Full name:</b>	H1.3 linker histone, cluster member
<b>Synonyms:</b>	H1D; H1.3; H1F3; H1s-2; HIST1H1D
<b>SwissProt:</b>	P16402
<b>IHC positive control:</b>	Human ovarian cancer and Human colorectal cancer
<b>IHC Recommend dilution:</b>	50-300



