

兔抗 MRPL28 多克隆抗体

中文名称：兔抗 MRPL28 多克隆抗体

英文名称： Anti-MRPL28 rabbit polyclonal antibody

别名： p15; MAAT1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MRPL28

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and some times in biochemical properties, which prevents easy reco

	gnition by sequence homology.
Applications:	ELISA, WB, IHC
Name of antibody:	MRPL28
Immunogen:	Fusion protein of human MRPL28
Full name:	mitochondrial ribosomal protein L28
Synonyms :	p15; MAAT1
SwissProt:	Q13084
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human thyroid cancer and human cervical cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	30 kDa
WB Positive control:	293T cells and human testis tissue, A375 cells and human normal kidney tissue
WB Recommended dilution:	1000-5000



