

IL23R 抗原（重组蛋白）

中文名称：IL23R 抗原（重组蛋白）

英文名称：IL23R Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 31-227 amino acids of human IL23R

技术规格

Full name:	Interleukin-23 receptor
Swissprot:	Q5VWK5
Gene Accession:	BC040720
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a subunit of the receptor for IL 23A/IL23. This protein pairs with the receptor molecule IL12RB1/IL12 Rbeta1, and both are required for IL23A signaling. This protein associates constitutively with Janus kinase 2 (JAK2), and also binds to transcription activator STAT3 in a ligand-dependent manner.