

兔抗 HOXC6 多克隆抗体

中文名称：兔抗 HOXC6 多克隆抗体

英文名称：Anti-HOXC6 rabbit polyclonal antibody

别名：CP25, HOX3, HOX3C, HHO.C8

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HOXC6

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the homeobox family, members of which encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, which are located on different chromosomes and consist of 9 to 11 genes arranged in tandem. This gene, HOXC6, is one of several HOXC genes located in a cluster on chromosome 12. Three genes, HOXC5, HOXC4 and HOXC6, share a 5' non-coding exon.
Applications:	ELISA, WB

Name of antibody:	HOXC6
Immunogen:	Synthetic peptide of human HOXC6
Full name:	homeobox C6
Synonyms :	CP25; HOX3; HOX3C; HHO.C8
SwissProt:	P09630
ELISA Recommended dilution :	2000-5000
WB Predicted band size:	27 kDa
WB Positive control:	Huvec cell
WB Recommended dilution:	500-2000

