

## 兔抗 MRPL3 多克隆抗体

中文名称：兔抗 MRPL3 多克隆抗体

英文名称： Anti-MRPL3 rabbit polyclonal antibody

别名： MRL3; RPML3; COXPD9

抗原： MRPL3

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein that belongs to the L3P ribosomal protein family. A pseudogene corresponding to this gene is found on chromosome 13q.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	MRPL3
<b>Immunogen:</b>	Synthesized peptide derived from Internal of human MRPL3.
<b>Full name:</b>	mitochondrial ribosomal protein L3
<b>Synonyms :</b>	MRL3; RPML3; COXPD9

<b>SwissProt:</b>	P09001
<b>WB Predicted band size:</b>	39 kDa
<b>WB Positive control:</b>	K562 cells lysate
<b>WB Recommended dilution:</b>	500-3000

