

兔抗 HOXC9 多克隆抗体

中文名称：兔抗 HOXC9 多克隆抗体

英文名称：Anti-HOXC9 rabbit polyclonal antibody

别名：homeobox C9; HOX3; HOX3B

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HOXC9

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, which are located on different chromosomes and consist of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXC genes located in a cluster on chromosome 12.
Applications:	ELISA, IHC
Name of antibody:	HOXC9

Immunogen:	Synthetic peptide of human HOXC9
Full name:	homeobox C9
Synonyms:	HOX3; HOX3B
SwissProt:	P31274
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer
IHC Recommend dilution:	100-200

