

兔抗 GATA3 多克隆抗体

中文名称：兔抗 GATA3 多克隆抗体

英文名称：Anti-GATA3 rabbit polyclonal antibody

级别：HDR; HDRS

抗原：GATA3

储存：冷冻 (-20℃)

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein which belongs to the GATA family of transcription factors. The protein contains two GATA-type zinc fingers and is an important regulator of T-cell development and plays an important role in endothelial cell biology. Defects in this gene are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia.
Applications:	ELISA, WB
Name of antibody:	GATA3
Immunogen:	Synthetic peptide of human GATA3
Full name:	GATA binding protein 3

Synonyms:	HDR; HDRS
SwissProt:	P23771
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	48 kDa
WB Positive control:	NIH/3T3 cell lysate
WB Recommended dilution:	500-2000

