

IRAK2 抗原（重组蛋白）

中文名称：IRAK2 抗原（重组蛋白）

英文名称：IRAK2 Antigen (Recombinant Protein)

别名：IRAK-2

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 111-210 amino acids of human IRAK2

技术规格

Full name:	interleukin 1 receptor associated kinase 2
Synonyms:	IRAK-2
Swissprot:	O43187
Gene Accession:	BC125184
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB.