

## 兔抗 GATA4(Phospho-Ser105) 多克隆抗体

中文名称：兔抗 GATA4(Phospho-Ser105) 多克隆抗体

英文名称： Anti-GATA4(Phospho-Ser105) rabbit polyclonal antibody

别名： TOF; ASD2; VSD1; TACHD

抗原： GATA4(Phospho-Ser105)

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a member of the GATA family of zinc-finger transcription factors. Members of this family recognize the GATA motif which is present in the promoters of many genes. This protein is thought to regulate genes involved in embryogenesis and in myocardial differentiation and function. Mutations in this gene have been associated with cardiac septal defects.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	GATA4(Phospho-Ser105)
<b>Immunogen:</b>	Peptide sequence around phosphorylation site of Serine 105 (P-V-S(p)-P-R) derived from Human GATA4.

<b>Full name:</b>	GATA binding protein 4
<b>Synonyms :</b>	TOF; ASD2; VSD1; TACHD
<b>SwissProt:</b>	P43694
<b>WB Predicted band size:</b>	45 kDa
<b>WB Positive control:</b>	293 cells lysates
<b>WB Recommended dilution:</b>	500-1000

