

IRF6 抗原（重组蛋白）

中文名称：IRF6 抗原（重组蛋白）

英文名称：IRF6 Antigen (Recombinant Protein)

别名：LPS; PIT; PPS; VWS; OFC6; PPS1; VWS1

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to C terminal 301 amino acids of human IRF6

技术规格

Full name:	interferon regulatory factor 6
Synonyms:	LPS; PIT; PPS; VWS; OFC6; PPS1; VWS1
Swissprot:	O14896
Gene Accession:	BC014852
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the interferon regulatory transcription factor (IRF) family. Family members share a highly-conserved N-terminal helix-turn-helix DNA-binding domain and a less conserved C-terminal protein-binding domain. The encoded protein may be a transcriptional activator. Mutations in this gene can cause van der Woude syndrome and popliteal pterygium syndrome. Mutations in this gene are also associated with non-syndromic orofacial cleft type

6. Alternate splicing results in multiple transcript variants.