

兔抗 HOXD3 多克隆抗体

中文名称：兔抗 HOXD3 多克隆抗体

英文名称：Anti-HOXD3 rabbit polyclonal antibody

别名：HOX4; HOX1D; HOX4A; Hox-4.1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HOXD3

反应种属：Human, Mouse

标记物：Unconjugate


克隆类型：Unconjugate

技术规格

Background:

This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, located on different chromosomes, consisting of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXD genes located at 2q31-2q37 chromosome regions. Deletions that removed the entire HOXD gene cluster or 5' end of this cluster have been associated with severe limb and genital abnormalities.

Applications:	ELISA, WB
Name of antibody:	HOXD3
Immunogen:	Synthetic peptide of human HOXD3
Full name:	homeobox D3
Synonyms :	HOX4; HOX1D; HOX4A; Hox-4.1
SwissProt:	P31249
ELISA Recommended dilution:	5000-10000
IHC Recommend dilution :	0
WB Predicted band size:	46 kDa
WB Positive control:	HEPG2 cell
WB Recommended dilution:	200-1000



mlbio 梅联生物
Good elisakit producers