

## IRS1 抗原（重组蛋白）

中文名称：IRS1 抗原（重组蛋白）

英文名称：IRS1 Antigen (Recombinant Protein)

别名：HIRS-1

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from N terminal 264 amino acids of human IRS1

技术规格

<b>Full name:</b>	insulin receptor substrate 1
<b>Synonyms:</b>	HIRS-1
<b>Swissprot:</b>	P35568
<b>Gene Accession:</b>	BC053895
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein which is phosphorylated by insulin receptor tyrosine kinase. Mutations in this gene are associated with type II diabetes and susceptibility to insulin resistance.