

兔抗 MRPL49 多克隆抗体

中文名称：兔抗 MRPL49 多克隆抗体

英文名称： Anti-MRPL49 rabbit polyclonal antibody

别名： NOF; NOF1; L49mt; C11orf4; MRP-L49

抗原： MRPL49

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. Pseudogenes corresponding to this gene are found on chromosomes 5q and 8p.
Applications:	WB
Name of antibody:	MRPL49
Immunogen:	Synthesized peptide derived from Internal of human MRPL49

Full name:	mitochondrial ribosomal protein L49
Synonyms :	NOF; NOF1; L49mt; C11orf4; MRP-L49
SwissProt:	Q13405
WB Predicted band size:	19 kDa
WB Positive control:	Hela cells and COLO205 lysates
WB Recommended dilution:	500-3000

