

兔抗 MRPL4 多克隆抗体

中文名称：兔抗 MRPL4 多克隆抗体

英文名称： Anti-MRPL4 rabbit polyclonal antibody

别名： L4mt; CGI-28

抗原： MRPL4

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. Sequence analysis identified alternatively spliced variants that encode different protein isoforms.
Applications:	WB
Name of antibody:	MRPL4
Immunogen:	Synthesized peptide derived from internal of human MRPL4.
Full name:	mitochondrial ribosomal protein L4
Synonyms :	L4mt; CGI-28

SwissProt:	Q9BYD3
WB Predicted band size:	35 kDa
WB Positive control:	COLO cells lysate
WB Recommended dilution :	500-3000

