

## ITGB6 抗原（重组蛋白）

中文名称：ITGB6 抗原（重组蛋白）

英文名称：ITGB6 Antigen (Recombinant Protein)

别名：AI1H

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 139-434 amino acids of human ITGB6

技术规格

<b>Full name:</b>	integrin, beta 6
<b>Synonyms:</b>	AI1H
<b>Swissprot:</b>	P18564
<b>Gene Accession:</b>	BC121178
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein that is a member of the integrin superfamily. Members of this family are adhesion receptors that function in signaling from the extracellular matrix to the cell. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. The encoded protein forms a dimer with an alpha v chain and this heterodimer can bind to ligands like fibronectin and transforming growth factor beta 1. Alternate s

licing results in multiple transcript variants.