

JUN 抗原（重组蛋白）

中文名称：JUN 抗原（重组蛋白）

英文名称：JUN Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to C terminal 250 amino acids of human JUN

技术规格

Full name:	jun proto-oncogene
Synonyms:	AP1; AP-1; c-Jun
Swissprot:	P05412
Gene Accession:	BC006175
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.