

## 兔抗 GBP5 多克隆抗体

中文名称：兔抗 GBP5 多克隆抗体

英文名称： Anti-GBP5 rabbit polyclonal antibody

别名： GBP-5

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： GBP5

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	GBP5 (guanylate binding protein 5), also known as GBP-T A antigen, is a 586 amino acid protein that localizes to the cytoplasmic side of the cell membrane. Belonging to the interferon (IFN)-inducible guanylate-binding protein (GBP) family, GBP5 may be involved in the inflammatory response or with cell proliferation and is suggested to have possible cancer-related functions. The gene encoding GBP5 maps to human chromosome 1, which spans 260 million base pairs, contains over 3,000 genes and comprises nearly 8% of the human genome.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	GBP5
<b>Immunogen:</b>	Fusion protein of human GBP5
<b>Full name:</b>	guanylate binding protein 5
<b>Synonyms :</b>	GBP-5
<b>SwissProt:</b>	Q96PP8
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human prostate cancer and human brain
<b>IHC Recommend dilution:</b>	50-200

