

KRT17 抗原（重组蛋白）

中文名称：KRT17 抗原（重组蛋白）

英文名称：KRT17 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 102-346 amino acids of human KRT17

技术规格

Full name:	keratin 17
Synonyms:	PC; K17; PC2; 39.1; CK-17; PCHC1
Swissprot:	Q04695
Gene Accession:	BC000159
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes the type I intermediate filament chain keratin 17, expressed in nail bed, hair follicle, sebaceous glands, and other epidermal appendages. Mutations in this gene lead to Jackson-Lawler type pachyonychia congenita and steatocystoma multiplex.