

兔抗 GCC2 多克隆抗体

中文名称：兔抗 GCC2 多克隆抗体

英文名称：Anti-GCC2 rabbit polyclonal antibody

抗 原：GCC2

别 名：GRIP and coiled-coil domain containing 2; REN53; GCC185; RANBP2L4

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a peripheral membrane protein localized to the trans-Golgi network. It is sensitive to brefeldin A. This encoded protein contains a GRIP domain which is thought to be used in targeting. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	GCC2
Immunogen:	Synthetic peptide of human GCC2
Full name:	GRIP and coiled-coil domain containing 2
Synonyms:	REN53; GCC185; RANBP2L4
SwissProt:	Q8IWJ2

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-300

