

KRT31 抗原（重组蛋白）

中文名称：KRT31 抗原（重组蛋白）

英文名称：KRT31 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 100-400 amino acids of human KRT31

技术规格

Full name:	Keratin-31
Synonyms:	HA1, Ha-1, hHa1, KRTHA1
Swissprot:	Q15323
Gene Accession:	BC114468
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the keratin gene family. As a type I hair keratin, it is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription.