

## LGALS2 抗原（重组蛋白）

中文名称：LGALS2 抗原（重组蛋白）

英文名称： LGALS2 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-132 amino acids of human LGALS2

技术规格

<b>Full name:</b>	galectin 2
<b>Synonyms:</b>	HL14
<b>Swissprot:</b>	P05162
<b>Gene Accession:</b>	BC059782
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a soluble beta-galactoside binding lectin. The encoded protein is found as a homodimer and can bind to lymphotoxin-alpha. A single nucleotide polymorphism in an intron of this gene can alter the transcriptional level of the protein, with a resultant increased risk of myocardial infarction.