

## 兔抗 HPRT1 多克隆抗体

中文名称：兔抗 HPRT1 多克隆抗体

英文名称：Anti-HPRT1 rabbit polyclonal antibody

别名：HPRT; HGPRT

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：HPRT1

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a transferase, which catalyzes conversion of hypoxanthine to inosine monophosphate and guanine to guanosine monophosphate via transfer of the 5-phosphoribosyl group from 5-phosphoribosyl 1-pyrophosphate. This enzyme plays a central role in the generation of purine nucleotides through the purine salvage pathway. Mutations in this gene result in Lesch-Nyhan syndrome or gout.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	HPRT1
<b>Immunogen:</b>	Fusion protein of human HPRT1
<b>Full name:</b>	hypoxanthine phosphoribosyltransferase 1

<b>Synonyms :</b>	HPRT; HGPRT
<b>SwissProt:</b>	P00492
<b>IHC positive control:</b>	Human prostate cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	25 kDa
<b>WB Positive control:</b>	Hela, 231, hepg2 and A549 cell
<b>WB Recommended dilution:</b>	500-2000



