

兔抗 MRPS34 多克隆抗体

中文名称：兔抗 MRPS34 多克隆抗体

英文名称： Anti-MRPS34 rabbit polyclonal antibody

别名： MRPS12; MRP-S12; MRP-S34

抗原： MRPS34

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology.
Applications:	WB
Name of antibody:	MRPS34

Immunogen:	Synthesized peptide derived from C-terminal of human MRPS34.
Full name:	mitochondrial ribosomal protein S34
Synonyms :	MRPS12; MRP-S12; MRP-S34
SwissProt:	P82930
WB Predicted band size:	26 kDa
WB Positive control:	Hela cells, Huvec cells, HT-29 cells, A549 cells and Jurkat cells lysates
WB Recommended dilution:	500-3000

