

兔抗 HABP4 多克隆抗体

中文名称：兔抗 HABP4 多克隆抗体

英文名称：Anti-HABP4 rabbit polyclonal antibody

别名：hyaluronan binding protein 4; IHABP4; IHABP-4; Ki-1/57; SERBP1L

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HABP4

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in nuclear functions such as the remodeling of chromatin and the regulation of transcription.
Applications:	ELISA, IHC
Name of antibody:	HABP4
Immunogen:	Synthetic peptide of human HABP4
Full name:	hyaluronan binding protein 4
Synonyms:	IHABP4; IHABP-4; Ki-1/57; SERBP1L
SwissProt:	Q5JVS0

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and Human liver cancer
IHC Recommend dilution:	40-200



