

RAF1 抗原（重组蛋白）

中文名称：RAF1 抗原（重组蛋白）

英文名称：RAF1 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 349-609 amino acids of human RAF1

技术规格

Full name:	Raf-1 proto-oncogene, serine/threonine kinase
Synonyms:	NS5; CRAF; Raf-1; c-Raf; CMD1NN
Swissprot:	P04049
Gene Accession:	BC018119
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is the cellular homolog of viral raf gene (v-raf). The encoded protein is a MAP kinase kinase kinase (MAP3K), which functions downstream of the Ras family of membrane associated GTPases to which it binds directly. Once activated, the cellular RAF 1 protein can phosphorylate to activate the dual specificity protein kinases MEK1 and MEK2, which in turn phosphorylate to activate the serine/threonine specific protein kinases, ERK1 and ERK2. Activated ERKs are pleiotropic effectors of cell physiology and play an important role in the control of gene expression involved in the cell division cycle, apoptosis, cell differentiation and cell mi

gration.