

## 兔抗 HP 多克隆抗体

中文名称：兔抗 HP 多克隆抗体

英文名称：Anti-HP rabbit polyclonal antibody

别名：BP; HPA1S; HP2ALPHA2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HP

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a preproprotein, which is processed to yield both alpha and beta chains, which subsequently combine as a tetramer to produce haptoglobin. Haptoglobin functions to bind free plasma hemoglobin, which allows degradative enzymes to gain access to the hemoglobin, while at the same time preventing loss of iron through the kidneys and protecting the kidneys from damage by hemoglobin. Mutations in this gene and/or its regulatory regions cause alpha hemoglobinemia or hypohaptoglobinemia.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	HP
<b>Immunogen:</b>	Synthetic peptide of human HP

<b>Full name:</b>	haptoglobin
<b>Synonyms :</b>	BP; HPA1S; HP2ALPHA2
<b>SwissProt:</b>	P00738
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control :</b>	Human thyroid cancer and Human gastric cancer
<b>IHC Recommended dilution:</b>	100-300

