

兔抗 MSH2 多克隆抗体

中文名称: 兔抗 MSH2 多克隆抗体

英文名称: Anti-MSH2 rabbit polyclonal antibody

别名: mutS homolog 2 FCC1; COCA1; HNPCC; LCFS2; HNPCC1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: MSH2

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Synthetic peptide of human MSH2
Applications:	ELISA, WB, IHC
Name of antibody:	MSH2
Immunogen:	Synthetic peptide of human MSH2
Full name:	mutS homolog 2
Synonyms:	FCC1; COCA1; HNPCC; LCFS2; HNPCC1
SwissProt:	P43246
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human colorectal cancer;
IHC Recommend dilution:	20-100

WB Predicted band size:	105 kDa
WB Positive control:	HEPG2 , 293T , SKOV3 , Hela and A172 cell lysates
WB Recommended dilution:	200-1000



