

## 兔抗 HRG 多克隆抗体

中文名称：兔抗 HRG 多克隆抗体

英文名称：Anti-HRG rabbit polyclonal antibody

别名：HPRG; HRGP; THPH11

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

抗原：HRG

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This histidine-rich glycoprotein contains two cystatin-like domains and is located in plasma and platelets. The physiological function has not been determined but it is known that the protein binds heme, dyes and divalent metal ions. It can inhibit rosette formation and interacts with heparin, thrombospondin and plasminogen. Two of the protein's effects, the inhibition of fibrinolysis and the reduction of inhibition of coagulation, indicate a potential prothrombotic effect. Mutations in this gene lead to thrombophilia due to abnormal histidine-rich glycoprotein levels.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	HRG

<b>Immunogen:</b>	Synthetic peptide of human HRG
<b>Full name:</b>	histidine-rich glycoprotein
<b>Synonyms :</b>	HPRG; HRGP; THPH11
<b>SwissProt:</b>	P04196
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human breast cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	50-200

