

## 兔抗 MSH4 多克隆抗体

中文名称：兔抗 MSH4 多克隆抗体

英文名称： Anti-MSH4 rabbit polyclonal antibody

抗 原： MSH4

储 存： 冷冻（-20℃）

宿 主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the DNA mismatch repair mutS family. This member is a meiosis-specific protein that is not involved in DNA mismatch correction, but is required for reciprocal recombination and proper segregation of homologous chromosomes at meiosis I. This protein and MSH5 form a heterodimer which binds uniquely to a Holliday Junction and its developmental progenitor, thus provoking ADP-ATP exchange, and stabilizing the interaction between parental chromosomes during meiosis double-stranded break repair.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MSH4
<b>Immunogen:</b>	Fusion protein of human MSH4

<b>Full name:</b>	mutS homolog 4
<b>SwissProt:</b>	O15457
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and human brain
<b>IHC Recommend dilution:</b>	50-200

